



Patient Instructions for receiving calls from iCode Connect™

Dear Patient,

Your RESmart® CPAP comes equipped with iCode®. iCode® is an advanced way for your compliance data to be communicated to your provider. Your insurance carrier requires compliance monitoring for your therapy to be covered. iCode Connect™ is the system that collects your iCode® through a simple phone call. Your provider will work with you on setting up your call schedule to retrieve your iCode®

**To access your iCode® on your RESmart®
follow these simple instructions.**

From the Stand-By Screen **(A)** press and release the Heater Button **(B)**.
(see pictures on below)

(NOTE there are 6 iCode® periods, iCode® 1 day, 7 day, 30 day, 60 day, 90 day and 182 day. Your provider will tell you which iCode® you will be requested to communicate)

The display will now read: **(C)**

Pressing the heater button **(B)** again, you will now view a code similar to screen **(D)**.

Each time you press the HEATER **(B)** button, the display will advance to the next iCode® period FOLLOWED by the associated iCode®.



3B Medical, Inc.

On the dates your provider has scheduled iCode Connect™ calls for you, you should access your iCode® from your RESmart® device using the instructions provided and write down the iCode® your provider is requesting from you on the morning the call is scheduled. The number of calls scheduled depends on your insurance carrier's requirements and or that of your physician. Your iCode® is a string of 16 numbers, you can either say your iCode® string or enter it using the keypad on your phone. The system is similar to reading your Credit Card Number over the phone during a telephone purchase.

First Scheduled iCode Connect™ call date is: ____ / ____ / ____ Calls is Scheduled for ____:____ AM/PM

iCode requested by your provider is circled 1 ; 7 ; 30 ; 60 ; 90 ; 182

____ _
____ _

Second Scheduled iCode Connect™ call date is: ____ / ____ / ____ Calls is Scheduled for ____:____ AM/PM

iCode requested by your provider is circled 1 ; 7 ; 30 ; 60 ; 90 ; 182

____ _
____ _

Third Scheduled iCode Connect™ call date is: ____ / ____ / ____ Calls is Scheduled for ____:____ AM/PM

iCode requested by your provider is circled 1 ; 7 ; 30 ; 60 ; 90 ; 182

____ _
____ _

Fourth Scheduled iCode Connect™ call date is: ____ / ____ / ____ Calls is Scheduled for ____:____ AM/PM

iCode requested by your provider is circled 1 ; 7 ; 30 ; 60 ; 90 ; 182

____ _
____ _

3B Medical, Inc.

21301 US Highway 27 | Lake Wales, FL 33859 | Phone: (863) 226-6285 | Fax: (863) 226-6284 | info@3Bproducts.com | www.3BProducts.com