






# RESmart® Product Comparison Chart

Manufacturer	3B-BMC	ResMed	ResMed	Respironics	Respironics
Model	RESmart® Auto	S9 Elite	S9 Escape	PR550	DS220S
Picture					
Pressure Range (cmH <sub>2</sub> O)	4-20	4-20	4-20	4-20	4-20
Exhalation Pressure Relief	Yes (RESlex)	Yes (EPR)	Yes (EPR)	Yes (A-Flex)	Yes (C-Flex)
Compliance Remote Access*	Yes	No	No	No	No
Delay Off**	Yes	No	No	No	No
Safe Low voltage Humidifier	Yes	Yes	Yes	Yes	Yes
Ramp (min)	0-60	0-45	0-45	0-45	0-45
Automatic On/Off	Yes	Yes	No	Yes	No
Leak Alert	Yes	Yes	No	Yes	Yes
Power Off Alert	Yes	No	No	No	No
Memory (day)	365+1	365	365	365+7	365
Sound Level (dB)	<30	<30	<30	<30	<30
AHI	Yes	Yes	No	Yes	No
Leak	Yes	Yes	No	Yes	No
Compliance	Yes	Yes	Yes	Yes	Yes

\* iCode® permits remote compliance data retrieval by phone. Report generation includes reporting of Apnea-Hypopnea Index, Mean Pressure, P95, Snore Index, Number Days Compliant (<4 Hours), Best 30 Day Adherence Score, and % High Leak Data.

\*\* Delay Off is a feature designed to protect the electronics of the device, as well as to dry a patient's air hose and interface.

3B Medical, Inc.  
 21301 US Highway 27 N • Lake Wales, FL 33859 • (863) 226-6285 • Fax (863) 226-6284 • info@3bproducts.com  
[www.3Bproducts.com](http://www.3Bproducts.com)