



3B iCodeConnect™

Innovative Solutions in Sleep Therapy™




Wireless Compliance Report for John Doe

Date generated: 06/08/2016

Patient information	Provider information	Physician information
Patient name: John Doe Patient ID: JDoe1960 Device serial: E2315123456 Gender: Male Date of birth: 1/01/1960 Therapy start: 10/20/2015 Street 1: Street 2: City: Winter Haven State: FL ZIP:	Organization: Test Street 1: 1234 Maple St. Street 2: City: Lake Wales State: FL ZIP: 33853 Phone:	Name: Organization: Street 1: Street 2: City: State: ZIP: Phone:

Compliance Summary

Start date	10/23/2015
End date	02/08/2016
Days in range	109
Days of use	101 / 109 days (92.7%)
Total usage in date range	18 days, 7 hours, 56 minutes
Average usage (on days used)	4 hours, 21 minutes
Average usage (on all days)	4 hours, 2 minutes
Maximum usage	9 hours, 33 minutes
Minimum usage	3 minutes
Days with usage >= 4 hours	59.6% (65 of 109 days)



Usage Summary

Mean pressure	5.9 cmH2O
Average P95	6.0 cmH2O
Highest daily average pressure	6.0 cmH2O
Lowest daily average pressure	4.0 cmH2O
Average high leak time	5 minutes
Average leak	45.9 LPM
Average AHI	2.5
Average CAI	1.1

Key

P95: device was at or under this pressure 95% of the time
 High leak time: amount of time spent blowing with a leak > 90 L/min.
 CAI: clear airway apnea index

iCodeConnect™, a cloud based, information management system, makes managing and treating your sleep disorder patients easier and more efficient. With iCodeConnect, you can manage patient records, develop patient compliance reports, track overall treatment progress of PAP therapy patients and much more.


iCodeConnect™ provides accurate timely data, when you want it, is HIPPA compliant and cloud based so your patient's information can be accessed by their physicians and other healthcare providers in real time, with no additional costs or delays.

Providing the best possible care takes a team effort. iCodeConnect™ allows for quick and easy sharing of vital patient information. The entire patient care team: respiratory therapists, physicians, clinicians and others have quick and easy access to the information they need, when they need it.

At 3B™ Medical, we understand the challenging demands of compliance reporting in an environment of declining reimbursements. For this reason, 3B™ offers more low or no-cost solutions for remote compliance data reporting than any other manufacturer, with options ranging from Cellular Modem, WiFi, QR Coding on smart phone apps, and Interactive Voice Recognition. At 3B™ we are leading the way, and innovating new ways to reduce cost and improve efficiency.

Whether it is capturing, managing, sharing or tracking your patient's data, 3B's iCodeConnect™ makes doing it easier, better and free!

Multiple Data Connectivity Options



WiFi

WiFi provides reliable, "real time" data from your Luna® devices at no cost. Luna® WiFi allows you and your referrals to capture and access your patient's detailed compliance data whenever it is needed. And because it is typically free of transmission fees, it is the modality of choice for today's cost challenges.

iCode® with IVR

iCode® with Interactive Voice Recognition is just one more way 3B™'s iCodeConnect™ is helping companies return profitability to their bottom line. With IVR, you have the peace of mind knowing that no matter what the challenges other forms of data capture may present, you can breathe easy knowing IVR is there to help.

Cellular

Cellular modem offers the ability to capture your patient's necessary compliance data from almost anywhere, all the time. 3B™ offers one of the lowest cost cellular options in the market and provides detailed data at your fingertips any time you need it.

QR Code

With iCode-QR, any smart device can simply "scan", capture and automatically upload patient sleep data with the use of the free smart phone app. QR and the app are always free and instantly upload your patient data. And because iCodeConnect™ is cloud based, you can access your patient's data anytime, anywhere.

SD Card

All 3B™ products come standard with SD Card technology. With iCodeConnect™, no special card reader, software or other proprietary tools need to be installed. Simply upload your patient's device SD Card into iCodeConnect™ and your patient's detailed report is at your fingertips, safely and securely.

Because iCodeConnect™ is cloud based, you can access your patient's data anytime, anywhere.

...or John, Doe

Provider information

Organization: Sample Co., Ltd.
Street 1: 1234 Maple St.
Street 2:
City: Lake Wales
State: FL
ZIP: 33853
Phone: (555) 555-5555

iCode Connect®

Filter by provider

Item: lake wales, FL

- Northside HC - Richmond, VA
- DEF Homecare - Small town, FL
- ABC Homecare - Anywhere, FL new store

Patient ID	Name	Sex	City	State	Clinician
1234567	Victor, Nanni	M	psl	FL	Greenbridge-Thomas, Desiree
1274285	Demo, Patient	M	psl	FL	Nanni, Carlos
1967122345678	jose, nanni	M	psl	FL	Nanni, Carlos
56473839390	DEMO, PATIENT413	M	psl	FL	Nanni, Carlos
573482384829	juan, nanni	M	psl	FL	Nanni, Carlos
7654789796	nanni, rena	M	psl	FL	Nanni, Carlos
8988756433	lucio, jimmy	M	psl	FL	Nanni, Carlos
90765090	lucio, jimmy	M	psl	FL	Nanni, Carlos
9999999	Mike, Demo	M	psl	FL	Nanni, Carlos
ADEMO1967	DEMO, ANNE	F	psl	FL	Nanni, Carlos

Showing 1 to 10 of 19 entries

Select date to view: October 5th, 2013

Usage data for October 5th, 2013