



## Technical Service Bulletin

<b>Date:</b> 03/15/2018	<b>Model(s):</b> Cirrus 5
<b>Subject:</b> High Altitude Adjustment of the Cirrus 5 Oxygen Concentrator	<b>TSB#:</b> CIR2018102

### High Altitude Adjustment of the Cirrus 5 Oxygen Concentrator

This TSB rescinds TSB No. CIR2017100 which instructed on a manual altitude adjustment at elevations above 4,500 sft. With the release of Cirrus 5 V. 4 (i.e. serial numbers following prefix MZJ5S-101251) the Cirrus 5 no longer requires an altitude adjustment. The current Cirrus 5 Oxygen Concentrator is designed to make this elevation adjustment automatically with no input needed from end user.

7

*\*This TSB does not take the place of thoroughly reading the User/Service Manuals.  
Please refer to the user manual/service manual for complete operation and maintenance information.*