

Cirrus 5



- CIR1005 Cirrus 5 Oxygen Concentrator
- CIR2530 Cirrus 5 Oxygen Concentrator Filter Pack
- CIR5550 Cirrus 6" Humidifier Extension Tube
- O2S2007 3B Standard 7' Oxygen Cannula
- O2U2007 3B Ultra Soft 7' Oxygen Cannula
- O2T2025 3B 25' Oxygen Tubing
- O2C2030 3B Oxygen Tubing Connector

The Cirrus 5 is one of the smallest and quietest stationary concentrators available in the market.

Oxygen purity sensor is standard with every unit.

It is easy to use, provides up to 5 LPM of continuous, therapeutic oxygen and is virtually maintenance free. It is designed to operate in even the most challenging environments and provide years of reliable continuous oxygen.

3B™ Medical, Inc.

203 Avenue A NW Suite 300 • Winter Haven, FL 33881
info@3bproducts.com • www.3bproducts.com

(863) 226-6285

Model	Cirrus 5
Operation voltage (V/Hz)	AC 120V/60Hz
Oxygen flow (L/min)	.5 LPM - 5 LPM
Oxygen concentration (%)	93% +/- 3%
Outlet pressure (PSI)	5 - 10
Alarm	Power failure, low & high pressure, temperature, low purity.
Sound level (dB(A))	< 40
LCD Display	Number of Sessions, System Pressure, Time per Session and Total Hours
O2 Purity Sensor	Digital
Electrical category:	Class II Type B
Net Weight (Lbs)	37
Dimensions:	14" (W) x 12" (D) x 21" (H)
(OCI) Oxygen Concentration Indicator	When Oxygen Purity is >/= 85% Blue lamp is lit. When oxygen purity is <85%, red lamp is on, indicating low purity

Innovative Solutions in Respiratory and Sleep