



Introducing Lumin[®] Bullet

The newest, most ingenious way of disinfecting your CPAP hose,
from the inside out with no chemicals

Powerful germicidal UVC light
destroys germs and pathogens in 60 seconds
stops formation of biofilm build up



Outer Dimensions:
215 x 64 x 64 mm
8.5" x 2.5" x 2.5"

Weight:
10.5 OZ.

Item No:
LM4000

Battery last up to 30 uses
on one charge,
USB charging cord included

Compatible for use
on most popular CPAP
tubing brands

The Lumin[®] Bullet
is not compatible with the
ResMed Climate Line heated hose.

 **3B[™] Medical, Inc.**

Add: 203 Avenue A NW Suite 300, Winter Haven FL 33881
Tel: 863-226-6285 Fax: 863-226-6284

3865_04.19_luminbullet

Info@3Bproducts.com
www.3Bproducts.com