


About The TravelPAP

The NEW Luna TravelPAP is a small and lightweight PAP device that offers you the freedom to get a good night's sleep no matter where you go!



The Luna TravelPAP Offers

- » Small and Lightweight size, 14 ounces
- » Light ring indicator gives peace of mind that the device is working correctly
- » Supply Reminders
- » LightTrip (US) Smartphone APP
- » Optional 12V DC Adapter for Car, Camper, Truck or RV!



Luna TravelPAP Specifications

Size	6.26 in x 2.60 in x 2.83 in
Weight	14.1 ounces
Modes	CPAP, AutoCPAP
Pressure Range	4-20 cmH2O (in 0.5cm increments)
Ramp	Auto, 0-60 mins
Sound Level	<30dBA*

Great For

- » Business Travelers
- » Leisure Travelers
- » Truck Drivers
- » Car and RV Camping
- » and More!

Part Numbers

LT2000	Luna TravelPAP Auto CPAP System
LT1020	Air Filter For Luna TravelPAP
LTC30	Luna TravelPAP Carrying Case
BT15MM	Breathing Tube 15 mm
LTP200	Luna TravelPAP AC Power Adapter
PTRLT01	Luna TravelPAP Filter Cap
LTPDC	Luna TravelPAP DC Power Adapter



www.reacthealth.com/sleep/pap-devices