



# REACT HEALTH CONNECT

TM

## New Cloud-Based Data Sharing Platform For Sleep and Respiratory Devices



Market-leading ability to drill into data, with extended search for physicians



Data from all sources consolidated in one location, including wireless, SD card, and QR+ code uploads



Enhanced graphic layout with intuitive, user-friendly navigation



Refreshed reporting, enhanced security, HIPAA compliant



Remote management of patient prescription




Care teams can remotely view trends, track compliance, and make informed treatment decisions

# Patient Data\*

Customizable patient data displays

- ✓ Configure date range
- ✓ See daily averages
- ✓ Confirm initial compliance
- ✓ Toggle between horizontal, vertical and details views




Gray Bar: No data available

Red Bar: Under 4 hours of usage


Green Bar: 4+ hours of usage

White Bar: Device turned on but not used that day

Click into any day to view color-coded activity graph



Click and drag over any section of graph to zoom in



# Compliance and Therapy Reports

Easily access specific patient records from intuitive patient menu


Patients standard view Extended search

Last name First name Patient ID Date of Birth Serial number  Filters 0 7 15

NAME ↑	BIRTH DATE ↓↑	THERAPY START ↓↑	DAYS ON SERVICE ↓↑	COMPLIANCE ↓↑	STATUS ↓↑
ALL, Sue	05/09/1971	05/27/2023	116	27%	Active
ben, ben	06/23/1990	10/20/2022	335	✓	Active
Bridget, Bridget	03/17/1971	03/22/2022	547	7%	Active
Czapla, Piotr	01/01/2001	01/27/2023	236	0%	Active
Daniel, Daniel	08/27/1965	02/06/2023	226	13%	Active
Kim, Kim	03/03/2022	06/01/2023	111	0%	Active
Mike, Mike	01/01/1970	09/04/2023	16	37%	Active
Patient, Patient	01/01/1970	12/01/2021	658	✓	Active
poulosteat, t	01/01/1970	08/23/2022	393	✓	Active

Showing 1 to 12 of 12 entries

Quick, easy PDF export of data rich, user-friendly reports



To learn more or begin your enrollment,  
scan the QR code below



React Health is an American sleep and respiratory device manufacturer based in Sarasota, Florida

## Respiratory



**V\*Home**  
Ventilator with  
High Flow  
Therapy

**VC+Pro**  
Ventilator  
with  
Integrated  
Cough-Assist



Remote Download Enabled  
with React DataLink accessory

## Masks



**Siesta**  
Full Face &  
Nasal Masks

**Rio II**  
Nasal Pillows  
Mask



## Sleep



**Luna G3 CPAP & APAP**



**Luna G3 BiLevel 25A**  
**Luna G3 BiLevel 30 VT S/T**

For additional information, contact your React Health representative or visit [reacthealth.com](https://reacthealth.com)



[Reacthealth.com](https://reacthealth.com)

© 2023 React Health, all rights reserved.

MI.MK.0382.REVA