

React Health PAP Mask Part Numbers:



Siesta Full Face Mask

SFF1000: Siesta Full Face Mask- All Size Kit
SFF1001: Siesta Full Face Mask- S
SFF1002: Siesta Full Face Mask- M
SFF1003- Siesta Full Face Mask- L

Siesta Nasal Mask



SNM1000: Siesta Nasal Mask- All Size Kit
SNM1001: Siesta Nasal Mask- S
SNM1002: Siesta Nasal Mask- M
SNM1003- Siesta Nasal Mask- L

Rio II Nasal Pillows Mask

RII1000: Rio II Nasal Pillows Mask- All Size Kit
RII1001: Rio II Nasal Pillows Mask- S
RII1002: Rio II Nasal Pillows Mask- M
RII1003- Rio II Nasal Pillows Mask- L



REACTHEALTH



Our Family Of CPAP Mask Interfaces



REACTHEALTH

React Health
PAP Masks available at:



MI.MK.0276.REV.B.

The Siesta

Full Face Mask

- No forehead brace allows for increased line of sight - even with your glasses on
- Self-adjusting ergonomic cushion design allows for excellent seal
- Padded support arms minimize irritation and marks around cheekbones and eyes
- Quick release clips allow simple attachment and removal of mask
- Easy change cushion cup snaps into place allows for ease of removal and cleaning
- Swivel elbow for active sleepers and allows for sleep in any position
- Air vent diffuser for a quiet nights sleep
- Anti-phyxia valve safety feature opens for fresh air in case of machine shut off



The Siesta

Nasal Mask

- Padded support arms minimize irritation and marks around cheekbones and eyes
- Quick release clips allow for simple attachment and removal of mask
- Easy change cushion cup snaps into place allows for ease of removal and cleaning
- Swivel elbow for active sleepers allows for any sleep position



The Rio II

Nasal Pillows Mask

- Lightweight, flexible and quiet
- Dual wall pillows with trampoline bounce for secure fit
- Split back strap allows for a better fit, while Velcro fasteners create added layer of adjustment
- Unique swivel for active sleepers

