



INTRODUCING
Stratus 5

All new 5lpm oxygen concentrator
from 3B™ Medical

- Slim line design for a less obtrusive profile
- 360° swivel wheels to effortlessly roll over all types of flooring, including carpet
- Quiet sieve bed exchange
- Designed for high reliability

STR1005 - Stratus 5 Oxygen Concentrator

 **3B™ Medical, Inc.**



Stratus 5



Easy Access

to internal components for servicing

Noise Suppression

Compressor is specially placed and padded to dampen noise levels

Simple

intuitive screen.



Model STR1005

Rated power W	320
Voltage	120V±12V,60HZ±1 HZ
Flow rate L/min	0.5~5
Concentration Rated flow	93%±5%
Sound pressure level dB A	≤45
Sound power level dB A	≤55
Net weight kg	15.5
Dimension in	15.3 × 9 × 23.6

FEATURES:

Abnormal temperature alarm; Low oxygen purity alarm; Power failure alarm; Timer; Pressure alarm; Low flow alarm; Thermistor fault alarm

When the nominal pressure of the oxygen output port is zero, the oxygen concentration is 93%±5%

Product Codes

Replacement Parts

STR1020 - Stratus 5 Replacement Hepa Paper Filter

STR1030 - Stratus 5 Replacement Sponge Filter

STR1040 - Stratus 5 Replacement Bacteria Filter



O₂ Accessories

O2S2007 - 3B Standard 7' Oxygen Cannula

O2U2007 - 3B Ultra Soft 7' Oxygen Cannula

O2T2025 - 3B 25' Oxygen Tubing

O2C2030 - 3B Oxygen Tubing Connector



3B Medical, Inc.

Info@3Bproducts.com www.3Bproducts.com