

POC COMPARISON CHART

MANUFACTURER

3B Medical

GCE Zen O

Inogen

Chart

Respironics



POC	Aer X	Life	G3	Freestyle 5	Simply Go Mini
Height	8.38"	9.25"	8.25"	10.7"	9.4"
Depth	2.64"	3.8"	3"	4.4"	3.6"
Length	7.22"	9.8"	8.75"	6.6"	8.3"
Weight (lbs) (including 1 Battery)	4.25	5.5	4.8	6.7	5
Hours of Use (1 battery)	Up to 4 hours on pulse setting of 2pm	Up to 4 hours on pulse setting of 2pm	Up to 4.7 hours on setting of 2pm	Up to 2.5 hours on a setting of 2pm	Up to 4.5 hours on setting of 2pm
DbA	40dBA at 21pm	37dBA at 21pm	39DbA at 21pm	48DbA at 21pm	43DbA at 21pm
Patient Changeable Sieve	Yes	Yes	Yes	No	No
Pulse Settings	1-5	1-5	1-5	1-5	1-5
Bolus	1100ml	1050ml	1050ml	1000ml	1000ml

