

POC COMPARISON CHART

MANUFACTURER

3B Medical

GCE

Inogen One®

Chart

Phillips

Precision



POC	Aer X	Zen-O lite™	G5	Freestyle® 5	SimplyGo Mini	Live Active Five®
Height	8.38"	9.25"	8.15"	10.7"	9.4"	8.4"
Depth	2.64"	3.8"	3.26"	4.4"	3.6"	3.2"
Length	7.22"	9.8"	7.19"	6.6"	8.3"	8.5"
Weight (lbs) (including 1 Battery)	4.25	5.5	4.7	6.7	5	5
Hours of Use (1 battery)	Up to 4 hours on pulse setting of 2	Up to 4 hours on pulse setting of 2	Up to 6.5 hours on setting of 2	Up to 2.5 hours on setting of 2	Up to 4.5 hours on setting of 2	Up to 6.5 hours on setting of 2
DbA	40DbA at 2	37DbA at 2	38DbA at 2	48DbA at 2	43DbA at 2	40DbA at 2
Patient Changeable Sieve	Yes	Yes	Yes	No	No	Yes
Pulse Settings	1-5	1-5	1-6	1-5	1-5	1-6
Bolus	1100ml	1050ml	1260ml	1000ml	1000ml	1000ml
Warranty	5 years	3 years or 15,000 hours	3 years	3 years	3 years	3 years

