

Why Replace?

Mask Cushions begin to wear out and lose shape, once this happens they tend to leak.

Mask Frames are usually made of plastic and small cracks can appear in areas of high stress, wearing out hooks, clasps and especially moving parts.

Mask Headgear loses elasticity and Velcro tabs begin to wear out and lose their "gripping" power.

Hoses can have small rips in the soft plastic leading to unknown leaks and loss of pressure, hoses

are also the hardest component to keep clean and dry and biofilm may develop.

Air Filters wear out and deteriorate, making them ineffective against particles or they can have blockages and not filter effectively.

Water Chamber plastic can be damaged by high heat or chemicals and minerals in water, causing the plastic to become discolored over time. Silicone seals and gaskets wear out making them less water tight and leaving space to harbor germs.

CPAP Wipes

| | |
|--------|-------------------------------------|
| WA100 | Fragrance-Free (Aloe) 80/pack |
| WC100 | (Citrus) 80/pack |
| M200BX | Travel-Fragrance-Free (Aloe) 30/box |
| M201BX | Travel (Citrus) 30/box |

CPAP Tubing

| | |
|----------|---|
| UN100 | Ultra Noir (black) ONLY SOLD BY CASE |
| SAP100 | Sapphire (bright blue) ONLY SOLD BY CASE |
| RUB100 | Ruby (red) ONLY SOLD BY CASE |
| AME100 | Amethyst (purple) ONLY SOLD BY CASE |
| TOP100 | Topaz (lumin blue) ONLY SOLD BY CASE |
| PEA100 | Pink Pearl (light pink) ONLY SOLD BY CASE |
| 3BCL1000 | 3B ComfortLine Heated Tubing Kit (w/power supply) |
| 3BCL1010 | 3B Replacement ComfortLine Heated Tubing (Hybernite Compatible) |

| Product Type | HCPCS for Reimbursement | Replacement Schedule |
|----------------------------------|-------------------------|----------------------|
| Full Face Mask | A7030 | 1 per 3 months |
| Full Face Mask Cushion | A7031 | 1 per month |
| Nasal Mask 00 | A7034 | 1 per 3 months |
| Nasal Mask Cushion | A7032 | 2 per 1 month |
| Nasal Pillow Interface | A7034 | 1 per 3 months |
| Nasal Pillow Replacement Cushion | A7033 | 2 per 1 month |
| Mask Headgear (all types) | A7035 | 1 per 6 months |
| Filter Disposable | A7038 | 2 per 1 month |
| Filter Non-Disposable | A7039 | 1 per 6 months |
| Tubing | A7037 | 1 per 3 months |
| Tubing Heated | A4604 | 1 per 3 months |
| Water Chamber | A7046 | 1 per 6 months |



DEDICATED

Affordable Quality Sleep,
Respiratory & Disinfection
Products



A Full Range of Sleep Products
for the Home



www.3Bproducts.com | Info@3Bproducts.com

www.iCodeConnect.com

MID-179/v1.0_4360_02.21_sleep

SLEEP PRODUCTS

Luna®

| | |
|----------|--|
| LG2A00 | Luna II AutoCPAP |
| LG2000 | Luna II CPAP |
| LQX2A00 | Luna II QX Auto CPAP (Internet Model) |
| LQX2000 | Luna II QX CPAP (Internet Model) |
| LG3700 | Luna G3 BPAP 25A |
| LG3800 | Luna G3 BPAP 30VT BiLevel ST w Target VT Function |
| WL300 | Luna II 4G LTE |
| WL100 | WiFi Module (Comp w Luna & Luna II) |
| LG1020 | Air Filter - White (Comp w Luna & Luna II) <i>ONLY SOLD BY CASE</i> |
| L8008BX | Luna II Reusable Filter - Black <i>ONLY SOLD BY CASE</i> |
| LPSU | Luna Power Supply (Comp w Luna & Luna II) |
| LG24-510 | Luna II Water Chamber |
| LC B30 | Luna Replacement Carry Bag |
| LG230 | Luna II Replacement Carry Bag |



Rio® II Nasal

| | |
|---------|---------------------------|
| RII1000 | All Size Kit (S,M,L) |
| RII3001 | Replacement Pillow-Small |
| RII3002 | Replacement Pillow-Medium |
| RII3003 | Replacement Pillow-Large |
| RII2001 | Replacement Headgear |

Siesta Nasal Mask

| | |
|---------|------------------------------|
| SNM1000 | All Size Kit (S,M,L) |
| SNM1001 | Small |
| SNM1002 | Medium |
| SNM1003 | Large |
| SNM3001 | Replacement Cushion - Small |
| SNM3002 | Replacement Cushion - Medium |
| SNM3003 | Replacement Cushion - Large |
| SNM2001 | Replacement Headgear |

Siesta Full Face Mask

| | |
|---------|------------------------------|
| SFF1000 | All Size Kit (S,M,L) |
| SFF1001 | Replacement Cushion - Small |
| SFF1002 | Replacement Cushion - Medium |
| SFF1003 | Replacement Cushion - Large |
| SFF3001 | Large |
| SFF3002 | Medium |
| SFF3003 | Small |
| SFF2001 | Replacement Headgear |



Viva Nasal Mask

| | |
|--------|------------------------------|
| Vi1000 | All Size Kit (S,M,L) |
| Vi1001 | Small |
| Vi1002 | Medium |
| Vi1003 | Large |
| Vi3001 | Replacement Cushion - Small |
| Vi3002 | Replacement Cushion - Medium |
| Vi3003 | Replacement Cushion - Large |
| Vi2001 | Replacement Headgear |

Numa Full Face Mask

| | |
|--------|------------------------------|
| NU1001 | Small |
| NU1002 | Medium |
| NU1003 | Large |
| NU3001 | Replacement Cushion - Small |
| NU3002 | Replacement Cushion - Medium |
| NU3003 | Replacement Cushion - Large |
| NU2001 | Replacement Headgear |

Elara Full Face Mask

| | |
|---------|------------------------------|
| ELA1001 | Small |
| ELA1002 | Medium |
| ELA1003 | Large |
| ELA6010 | Replacement Cushion - Small |
| ELA6020 | Replacement Cushion - Medium |
| ELA6030 | Replacement Cushion - Large |
| ELA2001 | Replacement Headgear |