

# STATIONARY CONCENTRATOR COMPARISON

MANUFACTURER

3B Medical

Caire

Philips

Invacare

DeVilbiss

AirSep

5 LITER  
(SIDE BY SIDE)



5L CONCENTRATORS	Stratus 5	Companion 5	EverFlo	Perfecto 2	DeVilbiss	VisionAire 5
MODEL	<b>STR1005</b>	15067013	1020000	IRC5P	525	AS098-1
PERFORMANCE						
Continuous Oxygen Flow (max LPM)	<b>5</b>	5	5	5	5	5
Oxygen Concentration	<b>87%-96%</b>	87%-95.5%	90%-96%	87%-95.6%	90%-96%	87%-95.5%
Oxygen Outlet Pressure (max psi)	<b>6 PSI - 10 PSI</b>	5 PSI	5.5 PSI	5 PSI	8.5 PSI	7 PSI
Operating Altitude (feet above sea level)	<b>N/A</b>	-1,253 to 9,878	0 to 7,500	0 to 8,000	N/A	N/A
SIZE						
Cubic Size (inches)	<b>3228</b>	3628	3277	4140	2016	3372
Width (inches)	<b>13.78</b>	12.5	15	15	14	14.1
Height (inches)	<b>21.26</b>	21.5	23	23	12	20.8
Depth (inches)	<b>11.02</b>	13.5	9.5	12	12	11.5
Weight (pounds)	<b>32</b>	36	31	43	12	30
NOISE LEVEL (dBA)	<b>≤40</b>	40	45	43	37	45
WARRANTY	<b>3 yrs.</b>	3 yrs.	3 yrs.	3 yrs.	3 yrs.	3 yrs.